

# SAFE-T-SPORT® WRAP AROUND HINGED KNEE STABILIZING BRACE

**Model No: 37-350**

**Sizes: XSmall - 3XLarge**  
**Color: Black**

## HOW TO SIZE

Measure four inches above  
center of kneecap.

SIZE	FITS
X-Small	12-13"
Small	14-15"
Medium	16-17"
Large	18-19"
X-Large	20-21"
XX-Large	22-23"
XXX-Large	24-25"

## INDICATIONS:

- Slight to moderate ligament sprains and tears
- Medial/lateral meniscal tears
- Arthritis
- Knee joint instability


## FEATURES & BENEFITS:

- **SPORTS NEOPRENE** with a soft, long loop lining for therapeutic warmth and comfortable wear.
- **WRAP AROUND DESIGN** is ideal for swollen or injured knees. Front wrap closures are quick and easy to apply.
- **METAL MEDIAL/LATERAL HINGES** in neoprene covered pockets allow for full flexion while providing stabilization and preventing hyperextension.
- **OPEN POPLITEAL** eliminates bunching and irritation behind the knee.
- **ADJUSTABLE SUPPORT STRAPS** for an improved fit that stays in place.



## APPLICATION INSTRUCTIONS:

- Open adjustable straps and hook and loop closures on brace.
- Position brace behind leg so that popliteal opening in back of the brace is centered on the crease behind the knee. For a proper fit, the top of the brace should be positioned approximately six inches above the center of the kneecap. If applied correctly, the medial and lateral hinges should be aligned with the prominent bones on the each side of the kneecap.
- Secure the hook and loop thigh closure, then the calf closure. The neoprene should provide a comfortable amount of compression and should not feel constrictive.
- Insert top strap through the top plastic ring and pull strap back over itself, secure with hook and loop closure. Repeat for calf strap.
- To remove brace, open straps and hook and loop closures. Remember to check the hinge alignment with the center of the knee each time the brace is re-applied.



## CARE INSTRUCTIONS:

Close straps. Hand wash in cold water with mild soap. Wrap in towel to remove excess water (do not twist or wring). Air dry away from heat. Do not tumble dry. Do not use bleach.

**WARNING:** Contact of neoprene material with skin may cause irritation. Persons with a known susceptibility to dermatitis or with broken skin should avoid contact with this material. If skin irritation develops from contact with this material, discontinue use immediately and consult a physician. Should not be worn by those with phlebitis or circulatory disorders. Not recommended for use during sleep.  
**NOTE:** Consult your physician if you experience pain, swelling or unusual reactions.

800-327-4110



www.flarthopedics.com

©2007 FLA ORTHOPEDICS, INC.