



AFO NIGHT WRAP

Model No: 58-600

Sizes: SM/MD, LG/XL

Color: Black

INDICATIONS:

- Plantar Fasciitis
- Achilles Tendonitis
- Arthritis
- Calf Muscle Cramping
- Heel and Arch Pain
- Post-Surgical Treatment of the Foot

Size	Men's Shoe	Women's Shoe
SM/MD	6 - 9	up to 8
LG/XL	9 ½ and up	8 ½ - 11 ½

US Patent D515,217

800-327-4110

The HealWell® AFO Night Wrap helps prevent nighttime contracture with countervailing tension to hold the foot neutral. This helps relieve pain and discomfort and speed healing in new or mild conditions.

PRODUCT FEATURES:

- Soft laminated foam material is breathable, lightweight, and comfortable for nighttime wear with a smooth Lycra lining.
- Slim, low profile means no bulk; compact wrap design is easy to sleep in.
- Non-slip pad for safe, occasional walking.
- Compression formed back holds the wrap in place and conforms around the ankle and heel for improved foot positioning.
- Dual tension straps provide countervailing tension to hold the foot in position and prevent plantar flexion. Straps criss-cross for less bulk.
- Tension straps work with a low profile, rigid stay on the plantar surface, functioning as a lever to maintain foot position.
- Fits left or right



HOW TO APPLY:

Unfasten the tension straps and open the wrap. Place the foot inside with the heel in the opening and the foot flat on the floor.




First, close the wrap over the top of the foot, followed by the ankle closure and finally the calf wrap closure. Ensure wrap fits snug and is not too loose or too tight.

The tension straps are used to hold the foot in position and provide a gentle stretch. Grab both ends of the tension straps, pull back to gently flex the foot. Cross the strap on the left around the top of the foot and attach to the back of the wrap on the right side. Cross the strap on the right over the first strap by attaching the loop patch over the small square hook. Attach the right strap to the back of the wrap on the left side. The straps should cross in an "X" on top of the foot.

Do not over-stretch. Reduce the tension if any pain or discomfort is experienced.

CARE INSTRUCTIONS:

Hand wash in warm water with mild soap. Air dry away from heat.



Website: www.flaorthopedics.com

©2004 FLA ORTHOPEDICS, INC.