



Graduated Compression Thigh Highs

15-20 mm Hg

Sheer Therapy® Hosiery

Lite Support Graduated Compression, 15-20 mm Hg

MODEL: H22XY

- Indications:
- Preventative relief
 - Relieves leg fatigue, ankle and foot swelling, and mild varicosities
 - Improves circulation

- Features:
- Sheer material
 - Reinforced toe
 - Lace top

Colors:



Informative
Color Packaging

Sold in Pairs



	SIZE	NUDE	BLACK
Thigh High	A	H2201	H2261
	B	H2202	H2262
	C	H2203	H2263
	D	H2204	H2264
	Q	H2205	H2265

See Sizing Chart on Page 60



20-30 mm Hg

Soft Fit Graduated Therapy

Moderate Support Graduated Compression, 20-30 mm Hg

MODELS: H38XY, H39XY

- Indications:
- Relief of chronic leg fatigue and heaviness; ankle and foot swelling; mild varicosities, edema, and post sclerotherapy

- Features:
- Microfiber softness
 - Moisture wicking
 - Elegant sheen
 - Feeling of tights with therapeutic benefits
 - Available in thigh high with lace top and mid thigh high with Uni-Band top

Colors:



Informative Color Packaging
Sold in Pairs



	SIZE	BARELY BEIGE	IVORY	BLACK
Thigh High w/Lace	Small	H3801	H3811	H3861
	Medium	H3802	H3812	H3862
	Large	H3803	H3813	H3863
	X-Large	H3804	H3814	H3864
	Queen	H3805	H3815	H3865
For Men & Women Mid Thigh w/Uni-Band	Small	H3901	—	H3961
	Medium	H3902	—	H3962
	Large	H3903	—	H3963
	X-Large	H3904	—	H3964
	XX-Large	H3905	—	H3965

See Sizing Chart on Page 60

ACTIVA® COMPRESSION HOSIERY & SOCKS



30-40 mm Hg

Surgical Weight Graduated Therapy

Heavy Graduated Compression, 30-40 mm Hg

MODEL: H42XY

- Indications:
- Relief of chronic leg fatigue and heaviness
 - Relief of ankle and foot swelling (edema); chronic venous insufficiencies; post sclerotherapy; post surgical; and prophylaxis thrombosis.

- Features:
- Uni-Band top
 - Closed toe

Color:



For Men & Women

Informative Color
Packaging

Sold in Pairs



	SIZE	BEIGE
Thigh High Closed Toe	Small	H4201
	Medium	H4202
	Large	H4203
	X-Large	H4204

See Sizing Chart
on Page 60

Graduated Compression Pantyhose 15-20 mm Hg

Sheer Therapy® Hosiery

Lite Support Graduated Compression, 15-20 mm Hg

MODEL: H21XY

- Indications:
- Preventative relief
 - Relieves leg fatigue, ankle and foot swelling, and mild varicosities
 - Improves circulation

- Features:
- Sheer material
 - Reinforced toe
 - Control top

Colors:



SIZE	WHITE	NUDE	BLACK
A	H2111	H2101	H2161
B	H2112	H2102	H2162
C	H2113	H2103	H2163
D	H2114	H2104	H2164
Q	—	H2105	H2165

See Sizing Chart on Page 60



15-20 mm Hg

Informative
Color
Packaging



20-30 mm Hg

Soft Fit Graduated Therapy

Moderate Support Graduated Compression, 20-30 mm Hg

MODEL: H37XY

- Indications:
- Relief of chronic leg fatigue and heaviness; ankle and foot swelling; mild varicosities, edema, and post sclerotherapy

- Features:
- Microfiber softness
 - Moisture wicking
 - Elegant sheen
 - Feeling of tights with therapeutic benefits
 - Control top

Colors:



20-30 mm Hg

Informative Color
Packaging

Sold in Pairs



SIZE	BARELY BEIGE	IVORY	BLACK
Small	H3701	H3711	H3761
Medium	H3702	H3712	H3762
Large	H3703	H3713	H3763
X-Large	H3704	H3714	H3764
Queen	H3705	H3715	H3765

See Sizing Chart
on Page 60

18 mm Hg

Anti-Embolism Stockings

Graduated Compression, 18 mm Hg

MODELS: H52XY-H54XY

- Indications:
- Thrombo Embolic Deterrent; intended for recumbent (non-ambulatory) patients as a means to help prevent venous stasis and embolism formation; and for pre- and post-operative procedures and mild varicosities without significant edema.

- Styles available:
- Knee High, Closed toe, White or Beige
 - Thigh High, Inspection Toe, "Bandless" Top, White
 - Thigh High, Closed toe, Uni-Band Top, White or Beige

Colors:



Informative Color
Packaging


Sold in Pairs



See Sizing Chart
on Page 60

	SIZE	BEIGE	WHITE
Thigh High Closed Toe	Small	H5201	H5411
	Medium	H5202	H5412
	Large	H5203	H5413
	X-Large	H5204	H5414
Thigh High Inspect. Toe	Small	—	H5211
	Medium	—	H5212
	Large	—	H5213
	X-Large	—	H5214
Knee High Closed Toe	Small	H5301	H5311
	Medium	H5302	H5312
	Large	H5303	H5313
	X-Large	H5304	H5314

FOLLOWING ARE GUIDELINES FOR HOW TO FIT AND MEASURE ACTIVA® PRODUCT



The products are designed to fit people who are mesomorphic in shape which typically accounts for 85% of the population.

Based on standard sizes and measurements, quick and easy to use weight and height charts have been developed for the Activa® products. These charts assume that based on height/weight combinations leg measurements are predictable and fit into specific ranges. If the fit is questionable, this guide has the exact measurements at the specific leg points needed to ensure the correct compression.

USE THIS MEASUREMENT POINT GUIDE

- g** - Circumference at top of thigh (about 2" below groin)
- f** - Circumference at center of thigh
- c** - Circumference at widest part of calf
- b** - Circumference at the narrowest part of the leg above the ankle

TABLE OF SIZE & COLOR

x	COLOR	y	SIZE
0	Beige / Nude	1	A / Small
1	Tan / Suntan	2	B / Medium
1	White/Ivory	3	C / Large
4	Navy	4	D / X-Large
5	Grey	5	Queen / XX-Large
6	Black		
7	Brown		

Note - Use these x and y values to indicate color and size for hosiery part numbers.

Sheer Therapy®
15-20 mm Hg • H21xy, H22xy, H23xy
MEASUREMENT IN INCHES

SIZE	ANKLE (b)	CALF (c)	MID-THIGH (f)	UPPER THIGH (g)
A	7 - 8"	11 - 12 ½"	16 ¼ - 18 ¼"	17 ¼ - 19 ½"
B	8 - 9 ½"	12 ½ - 14 ¼"	18 ¼ - 19 ¾"	19 ½ - 21 ¾"
C	9 ½ - 11"	14 ¼ - 16"	19 ¾ - 21 ¾"	21 ¾ - 24 ½"
D	11 - 12 ½"	16 - 17 ¾"	21 ½ - 24 ¾"	24 ½ - 27 ½"
Q	12 ½ - 14"	17 ¾ - 19 ¾"	24 ½ - 28"	27 ½ - 31"

Anti-Embolism Stockings
18 mm Hg
Hospital Style, Thigh High, Inspection Toe
Model H521y
CIRCUMFERENCE IN INCHES

SIZE	CALF. CIRCUM. (c)	LENGTH
Small	10 - 12"	28" or less
Medium	11 - 13"	28 - 30"
Large	12 - 15"	30 - 32"
X-Large	13 - 17"	30 - 34"

Soft Fit Graduated Therapy
20-30 mm Hg • H37xy, H38xy, H39xy, H30xy
MEASUREMENT IN INCHES

SIZE	ANKLE (b)	CALF (c)	THIGH (f)	WAIST	LENGTH
Small	6 ¼ - 7 ½"	10 - 13"	15 - 18"	Up to 30"	21 - 27"
Medium	7 - 8 ½"	11 ½ - 14"	17 ½ - 20 ½"	Up to 33"	23 - 29"
Large	8 ½ - 10"	12 ½ - 15"	20 - 23"	Up to 36"	25 - 31"
X-Large	10 - 11 ½"	13 ½ - 16"	22 ½ - 25 ½"	Up to 38"	27 - 33"
Q/X-Large	11 - 12 ½"	14 ½ - 17"	25 - 28"	Up to 44"	29 - 35"

Hosiery Style, Closed Toe
Model H53xy, H520y, H541y
CIRCUMFERENCE IN INCHES

SIZE	ANKLE (b)	CALF (c)
Small	6 ½ - 8"	10 ½ - 14 ½"
Medium	8 - 9 ½"	11 ½ - 15 ½"
Large	9 ½ - 11"	12 ½ - 17"
X-Large	11 - 12 ½"	13 ½ - 17 ½"

Surgical Weight Graduated Therapy
30-40 mm Hg • H42xy, H43xy, H44xy
MEASUREMENT IN INCHES

SIZE	ANKLE (b)	CALF (c)	MID-THIGH (f)	UPPER THIGH (g)
Small	6 ¼ - 7 ¼"	11 - 13 ½"	15 - 19"	16 - 20 ½"
Medium	7 ½ - 9"	13 - 15 ½"	18 ½ - 22"	20 - 24 ½"
Large	9 ¼ - 10 ¾"	15 - 17 ½"	21 ½ - 25"	24 - 29 ½"
X-Large	11 - 12 ½"	17 - 19 ½"	24 - 30"	29 - 33 ½"

SIZE	Men's & Women's Dress Socks		Assuming the following ankle and calf measurements in inches		Athletic Support Socks	
	MEN'S SHOE	WOMEN'S SHOE	ANKLE (b)	CALF (c)	MEN'S SHOE	WOMEN'S SHOE
Small	7 or Smaller	5 or Smaller	6 ¼ - 7 ¼"	11 - 13 ½"	7 or Smaller	7 - 9
Medium	7 ½ - 10	5 ½ - 7 ½	7 ½ - 9"	13 - 15 ½"	7 ½ - 10	9 - 12
Large	10 ½ - 12	8 - 10 ½	9 ¼ - 10 ¾"	15 - 17 ½"	10 ½ - 12	—
X-Large	12 ½ or Larger	11 or Larger	11 - 12 ½"	17 - 19 ½"	12 ½ or Larger	—

SOCK SIZING CHART

ACTIVA® COMPRESSION HOSIERY & SOCKS

©2004 FLA ORTHOPEDICS, INC., MIRAMAR, FL 33025