

IMPROVE
PATIENT COMPLIANCE!



PLANTAR FASCIITIS
NIGHT SPLINT
Resting Comfort Slipper™

MODEL: 58-500/501

The HealWell® Cub™ Resting Comfort Slipper™ is designed to comfortably position the foot in the controlled amount of dorsiflexion to provide a gentle stretch across the plantar fascia, Achilles tendon and calf muscles. This gentle stretch helps reduce the muscle and ligamentous contracture, inflammation and associated pain - clinically proven to reduce symptoms of plantar fasciitis. This action also helps relieve many other conditions of the lower extremities including foot drop, Achilles tendinitis, and painful night calf muscle cramping. Color: Black. Specify right or left.



US Patent D497,207; D495,058
US & Foreign Patents Pending

- Low profile, lightweight with less bulk for superior comfort.
- Open toe "slipper" is comfortable against the foot and easy to sleep in. Cool and comfortable to wear.
- Lateral side support with opening ensures no pressure on the sensitive heel, Achilles tendon, or bony ankle malleolus.
- Adjustable hinge can be positioned from 0 to 10 degrees depending on the amount of stretch required or desired.
- Built in toe wedge for more effective stretch across the plantar fascia.
- Ankle strap is perfectly positioned to hold the foot and heel in place for a more effective stretch.
- Easy and quick to apply.
- Cool, breathable foam material.



INFORMATIVE INSERT - Packaged with each product detailing information on plantar fasciitis and complete care & use instructions.

UNIQUE, ATTRACTIVE PACKAGING - Triangular shaped box with high impact full color graphic requires minimal shelf space.



SIZE	RIGHT EA ORDER #	LEFT EA ORDER #	CASE (4/CS) ORDER #	Fits Women's Shoe	Fits Men's Shoe
Small	58-5004	58-5014	HC5850SM	7½ and under	6½ and under
Medium	58-5005	58-5015	HC5851MD	8 - 11½	7 - 10½
Large	58-5006	58-5016	HC5852LG	12 and up	11 - 13

*4 Unit Case Contains 2 Right & 2 Left Per Size



ADJUSTABLE DORSIFLEXION

Adjustable hinge can be positioned from 0 to 10 degrees, depending on the amount of stretch required or desired.

Preset at most common 5 degree setting.

Indicated to Treat:
 Plantar Fasciitis • Achilles Tendinitis • Heel Spurs
 Muscle Tightening • Runner's Cramp



U.S. Patent
 D434,504



Our Standard Plantar Fasciitis Night Splint

HealWell® Plantar Fasciitis Night Splint

MODEL: 58-409

Comfortably positions the foot in a controlled amount of dorsiflexion. Adjustable bilateral dorsiflexion straps enable the application of prolonged stretch to the soft tissues of the Plantar Fascia and the Achilles tendon. Both the heel orientation and tibial straps are easy to adjust and use with push lock buckles. Padded straps for increased comfort. Made of high density plastic with ventilated calf area for patient comfort and breathability. Removable, soft liner. Removable foam wedge to provide additional toe extension. Comfortable C3™ Tri-Laminate material.

Universal design fits both left and right.



C3™ Tri-Laminate Liner
 for Cool and
 Comfortable Wear



ORDER #	SIZE	MEN'S SHOE SIZE	WOMEN'S SHOE SIZE
58-409400	Small	Under 6	5 - 7
58-409500	Medium	6 - 8	7½ - 9½
58-409600	Large	8½ - 10½	10 - 12
58-409700	X-Large	11 - 13	Over 12
58-409800	XX-Large	Over 13	—



Ankle Foot Orthosis/Foot Drop Splint

MODEL: 58-320

Indicated to correct foot drop during ambulation in patients with neurological deficit. Made of durable and resilient polypropylene. Fits into standard shoe. Can be easily trimmed or molded with heat. Forms at 310°-325°F. Adjustable calf band with padding and adjustable closure. Specify right or left.

RIGHT ORDER#	LEFT ORDER#	SIZE	MEN'S SHOE	WOMEN'S SHOE
58-320412	58-320421	Small	3 - 6	5 - 8
58-320519	58-320528	Medium	6 - 9	8 - 11
58-320616	58-320625	Large	9 - 11+	—