

# THERALL™ JOINT WARMING KNEE SUPPORT

**Model No: 53-702**

**Sizes: Small - XX-Large**  
**Color: Beige**

This elastic knee support provides light compression for comfortable support. Ideal for arthritis, weak or aching knees, overuse or repetitive use injuries and mild strains.

- Four-way stretch for light compression - reduces swelling
- Promotes healing by increasing circulation around the tender joint
- Improves knee mobility - allows for faster return to daily activities
- Use with Therall™ Therapeutic Heat Patch (sold separately) for over 12 hours of soothing heat therapy.



**Wearing Instructions:** Slide leg into support and pull up with pocket centered over knee.

**Care Instructions:** Hand wash in mild water with gentle soap. Air dry. Do not bleach or tumble dry. Ensure completely dry before reapplying.

**Material Content:** 12% Nylon, 25% Polyester, 33% Spandex, 30% Cotton, **LATEX FREE.**

**SIZING:** Measure around center of knee. Fits right or left.

<b>Size</b>	<b>Fits</b>
Small	11 - 13"
Medium	13 - 15"
Large	15 - 17"
X-Large	17 - 19 ½"
XX-Large	19 ½ - 22"

*Contains unique ceramic fibers that retain heat and slowly reflect it back into the joint and surrounding tissues. The result is soothing relief and therapeutic heat penetrated deep into the aching knee, muscles and surrounding tissues.*

**800-327-4110**



Website: **[www.flaorthopedics.com](http://www.flaorthopedics.com)**

©2008 FLA ORTHOPEDICS, INC.