



TRAINER ATHLETIC SHOES

For Men:

Model No: R10

Sizes: 7 - 12, 13, 14
(12 ½, 13 ½ not available)

Widths: M, W, XW

Color: White

For Women:

Model No: R50

Sizes: 6 - 10, 11
(10 ½ not available)

Widths: M, W, XW

Color: White


BENEFITS:

Properly fitted therapeutic shoes are critical to protecting sensitive feet. Shoes designed specifically for diabetic feet help prevent blisters, sores and skin breakdown.

Certain features not found in ordinary shoes, such as extra depth and width, help provide a proper fit and prevent the development of unwanted foot complications.

PRODUCT FEATURES:

- Breathable, flexible full grain leather upper
- Completely seamless interior with soft padding and breathable pique fabric interior
- Lining fabric upper wicks away moisture and minimizes heat build-up
- Extra depth, wide toe box and contoured last
- Padded heel collar reduces heel slippage and provides comfort
- Easy one-hand closure with traditional lace appearance
- Removable insole to accommodate custom orthotics
- Includes premium tri-layer removable insoles – two durometer polyurethane mold retains compression set and adds cushioning, shock blocker heel cushion
- EVA mid-sole filled with millions of nitrogen micro bubbles provides cushioning, comfort and extra shock absorption
- Lightweight, non-skid rubber sole for flexibility and non-slip traction



CARE INSTRUCTIONS:

Clean shoes with leather cleaner designed for soft leather. Shine with soft cotton cloth.

USE INSTRUCTIONS:

Wear shoes and insoles as prescribed by your medical professional. Shoes should be worn at all times, even indoors, to keep feet protected. Rotate insoles on a regular basis (every 3-4 months) to ensure proper pressure relief. *Relano Therapeutic Shoes are only one element in the overall diabetes management program. If any problems occur (discomfort, swelling, blisters, etc.) return immediately to your supplier to have the fitting personnel check for proper fit. If you notice a wound on the foot that does not begin to heal within one day, contact your medical professional immediately.*

FOOT CARE TIPS:

- Always keep your feet protected from the elements by wearing socks and shoes, even indoors.
- Wear therapeutic/diabetic shoes specially designed for sensitive, at-risk feet.
- Wear padded, seamless socks, such as PressureLite® Diabetic Socks, with shoes to help protect the feet
- Wash, but **do not soak** your feet in warm water each day and completely dry your feet, especially in between the toes.
- Check your feet for blisters, sores, cuts and swelling every day. Contact your healthcare provider if you notice a wound that does not begin to heal within one day.
- Never try to remove corns or calluses yourself. Speak with your healthcare provider about treatment.

800-327-4110



Website: www.flaorthopedics.com

©2006 FLA ORTHOPEDICS, INC.