


SOFT POINT® 3/4 FLEXIBLE CORK ORTHOTICS

Model No: F64

Sizes: Small - XLarge

Color: Blue

These 3/4 orthotics provide superior support to the arch of the foot and relieve the pain associated with many foot conditions. Deep contours cradle the heel and create full arch contact to comfortably stabilize the foot. The extremely lightweight and durable Cork/EVA composite material molds to the shape of the foot with the body's natural heat and weight, creating a customized feel.



Ideal for: Relief from Plantar Fasciitis; Forefoot, Arch and Heel Pain; Back and Joint Ache; Weak or Fallen Arches; Aching and Fatigued Feet; and Prevention of Pronation, or the inward rolling of the ankles.

FEATURES:

- **Extremely thin, flexible, lightweight base** is super lightweight, durable and shock absorbing.
- **Cork/EVA material molds to natural contours of the foot** to give a comfortable, custom-like fit and absorb shock.
- **Biomechanically designed** to provide longitudinal support across the foot and full support to the arch.
- **Breathable fabric** reduces perspiration, allowing for extended, comfortable wear and reduces shear to minimize blisters or calluses.
- **Deep heel cup** cradles the heel and helps stabilize the foot.
- **High side rims** allow for full contact with the arch, stabilizing the feet and preventing ankles from rolling inward.
- **3/4 length and moderate width** relieves pressure on the forefoot and accommodates a broad range of shoes – from walking and athletic shoes to casual and dress shoes. Can easily be moved from shoe to shoe.

HOW TO APPLY:

Place in shoe and trim edges to fit, if necessary. Use in closed toe shoes only. Always wear with socks or stockings. Can be easily transferred from one shoe to another.

CARE INSTRUCTIONS:

Wash in warm water with mild soap as needed. Air dry away from heat.

SIZING: Size according to US shoe size.

Size	Men's Shoe	Women's Shoe
Small	—	5 - 7
Medium	6 - 7	8 - 10
Large	8 - 10	—
X-Large	11 - 13	—

800-327-4110



Website: www.flaorthopedics.com

©2004 FLA ORTHOPEDICS, INC.