

# STANDARD NEOPRENE KNEE SLEEVE

Model No:

37-373  
Open Patella

37-374  
Closed Patella

Sizes: XSmall - XLarge  
Color: Navy

## INDICATIONS:

- Arthritis
- Meniscal Tears
- Joint Effusion
- Patellar Instability
- Knee Contusions
- Quad Weakness
- Knee Sprains/Strains

## BENEFITS:

These knee sleeves not only help to limit swelling, but also protect the injured knee. The neoprene generates heat, which provides increased warmth to the knee, thereby increasing circulation and decreasing pain. The thickness of neoprene also helps restrict excess motion and provides added support.

## PRODUCT FEATURES:

- Four way stretch for proper fit and comfort
- Sports neoprene for therapeutic warmth
- Easy slip-on style
- Open or closed patella

## APPLICATION:

Insert foot into knee support and slide the support up and over the knee. For proper fit, the top of the support should be positioned approximately four inches above the center of the kneecap.

## CARE INSTRUCTIONS:

Hand wash in cold water with mild soap before wearing. Wrap in towel to remove excess water (do not twist or wring). Air dry away from heat. Do not tumble dry. Do not bleach.

## SIZING:

Measure 4" above the center of the kneecap.

Size	Fits Knees
XSmall	12 - 13"
Small	14 - 15"
Medium	16 - 17"
Large	18 - 19"
XLarge	20 - 21"

## WARNING:

*Contact of this material with skin may cause irritation. Persons with a known susceptibility to dermatitis or with broken skin should avoid contact with this material. Anyone who develops skin irritation from contact with this material should discontinue its use immediately and consult a physician. Should not be worn by persons with phlebitis or circulatory disorders. Not recommended to be worn while sleeping.*

