

NEOPRENE PULL-ON THUMB SUPPORT

Model No: 25-130

Sizes: Small - X-Large

Color: Black or Beige

INDICATIONS:

- Arthritis
- Bursitis
- Tendonitis
- Repetitive motion injuries
- Soft tissue injuries
- General muscle soreness

BENEFITS:

This neoprene thumb support provides warmth and uniform compression to aching or injured thumbs and wrists. Compression helps reduce pain, inflammation, and swelling while providing support to soft tissues. Comfortable design allows for maximum hand use.



COMFORT BENEFITS:

- Full Thumb Coverage – Neoprene applies soothing warmth to all joints of the thumb.
- Contoured Shape – allows for full use of the hand.
- Soft Nylon Terry Lining – Works with soft stitching to give a seamless feel.

WEARING INSTRUCTIONS:

Simply pull thumb support on like a glove with the thumb through the smaller opening. The support is reversible to fit either right or left. Support should fit snugly but not so tight it digs into skin.

CARE INSTRUCTIONS:

Hand wash in cold water with mild soap. Wrap in a towel to remove excess water (do not twist or wring). Air dry away from heat. Do not tumble dry. Do not bleach.

SIZING: Measure circumference of wrist. Fits right or left.

<i>Size</i>	<i>Fits Wrists</i>
Small	5 ½ - 5"
Medium	6 - 6 ½"
Large	6 ½ - 7"
X-Large	7 - 8"

800-327-4110



Website: **www.flaorthopedics.com**

©2004 FLA ORTHOPEDICS, INC.